

ChemNeo

Semi-auto Chemistry Analyzer



General Features

- 7 inches touch screen display
- Long lifetime halogen lamp up to 2000 hours
- Permanent quartz flow cell
- Reaction curve real time display
- Thermal printer, Auto & Manually printing setting
- Flow cell & Cuvette modes
- Build in Incubator

Specification

- Principle: Colorimetry, Turbidimetry
- Wavelength: 340nm,405nm,450nm,505nm,546nm, 578nm,630nm and one free positions are available
- Light source: Halogen lamp,6V 10W
- Built-in Incubator : 8 + 2 Positions
- Display:7 inches touch screen
- Aspiration volume:100 ~ 5000uL
- Test method: End point, Fixed time, Kinetics, SUB, DIV, ABS
- Reagent: Open reagent system
- Absorbance variation ≤ 0.005 Abs
- Absorbance range: 0 ~ 3.000Abs
- Carry over $\leq 1\%$

Test Items

Liver	ALT, AST, ALP, GGT, TP, Alb, TBIL, DBIL
Myocardial enzyme	LDH, CK, CK-MB, α -HBDH
Lipids	CHO, TG, HDL_C, LDL_C
Lipoprotein	Apo A1, Apo B, Lp(α)
Kidney	BUN, CRE, UA
Others	CO ₂ -CP, AMS, LIP, GLU, HbA1C, IgA, IgG, IgM, C3, C4, CRP etc.

HETO Shenzhen Heto Medical Tech Co.,Ltd.

Address: 4/F Building No.4, Shenzhen Environmental Protection Industrial Park,
No.31 Makan South Road, Nanshan District, Shenzhen, 518055, P.R.China.

Tel: +86-755-28935967 / Fax: +86-755-28937587

Email: sales@hetomedical.com/ Web: www.hetomedical.com